Restorative Torque Indicator

A Torque Indication Wrench
for the application of 20 and 35Ncm
of tightening torque

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es

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pt

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BIOMET *3i* Restorative Torque Indicator

A Torque Indication Wrench for the application of 20 and 35Ncm of tightening torque

IMPORTANT PRODUCT INFORMATION

Instructions for Use

AVAILABILITY

The BIOMET *3i* Restorative Torque Indicator is packaged individually, in a non-sterile nylon pouch. Sterilization is possible by standard steam autoclaving. Please refer to the sterilization instructions section.

INDICATIONS FOR USE

The Restorative Indicator is indicated for 20 or 35Ncm of applied torque for the appropriate tightening of various posts and screw types that secure an abutment to a dental implant.

CONTRAINDICATIONS

There are currently no known contraindications for the use of the Restorative Torque Indicator when used for the indicated purpose of tightening dental abutment posts and abutment screws to the recommended torque of 20 or 35Ncm.

PRODUCT DESCRIPTION

The BIOMET *3i* Restorative Torque Indicator is a wrench with an infinite position, friction clutch drive mechanism that allows for the application of 20 and 35Ncm of tightening torque to dental abutments posts and abutment screws. The wrench is designed to be used with various interchangeable driver tips with the standard dental latch type design.

INSTRUCTIONS

Before Starting:

Before using the Restorative Torque Indicator, it is **important to note** that this wrench is only for the final application of 20 or 35Ncm of tightening torque. The selected screw or post should be lightly hand tightened, and proper seating should be confirmed radiographically, before the application of the final tightening torque with the wrench.

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


Securing the driver tip into the Restorative Torque Indicator:

Fig 1.
Select the appropriate dental latch type driver tip and secure it into the driver.

To deliver the desired torque:



Having selected and secured the appropriate driver tip into the Restorative Torque Indicator wrench, apply tightening torque to the screw or post in the normal clockwise direction. The wrench's friction clutching system allows the wrench to drive in the clockwise direction and slip in the counterclockwise direction.

Important Note: The Restorative Torque Indicator is for the indication of 20 and 35Ncm of applied torque only. Given the infinite clutch drive mechanism, it is possible to deliver torque in excess of 35Ncm. The application of loading forces beyond 35Ncm of torque can result in permanent damage to the unit and will void the warranty.

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

PRECAUTIONS

Failure to properly seat the driver tip fully into the wrench could result in premature tip disengagement. If this should occur intraorally, inadvertent aspiration or ingestion of the driver tip by the patient is possible.

United States federal law restricts the sale, distribution or use of this device to, by, or on the order of a licensed dentist or physician.

CAUTIONS

The Restorative Torque Indicator wrench is calibrated at time of manufacture. It is critical to the continued accuracy of the device that it only be used for the indication of 20 or 35Ncm of applied torque for the tightening of various post and screw types that secure an abutment to a dental implant. Any misuse of the device such as over torquing, or abuse such as excessive shock from dropping, can result in damage and incorrect torque indication readings. It is the ultimate responsibility of the Physician to insure that the device is properly maintained.

MAINTENANCE

Cleaning:

Scrubbing the Restorative Torque Indicator wrench with a soft plastic brush to remove visible contamination after each use is the only routine maintenance that is recommended or possible.

Lubrication:

None needed. Do not attempt to lubricate device.

Sterilization Instructions:

1. Remove contents from packaging prior to sterilization. Product is sold non-sterile. Sterilize prior to first use.
2. Re-sterilize after each patient exposure.
3. Steam autoclaving is the recommended means for processing this device. Effectiveness of any autoclave equipment or process is directly dependent upon numerous factors beyond the Company's control including: wrapping techniques employed; brand, model and condition of autoclave equipment; care and maintenance techniques employed; and operator knowledge and experience.

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4. Before autoclaving: **Ensure** all sterilization equipment is in proper working order as specified by the manufacturer; **Ensure** all equipment manufacturer's instructions are properly followed by trained and qualified personnel; **Ensure** actual cycle employed has been properly validated by the user for the device(s)/load configuration being processed; and, use of appropriate sterilization indicator devices for each autoclave cycle.

The Company suggests as a starting point for validation, the following process conditions:			
Materials	Temperature	Cold Start	Hot Start
Unwrapped instruments	134°C	25 min.	15 min.
Wrapped instruments	134°C	25 min.	20 min.
Instrument packs	134°C	35 min.	25 min.

*The Company cannot anticipate all possible equipment processes and/or conditions that may be encountered. The above suggested operation conditions are to be considered as a starting point for determination of the overall process capability, without regard for type or condition of equipment used or methods, techniques or practices employed by the user. Use of proper sterilization indicator devices is strongly recommended.

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